

دائرة اللوازم والمشتريات

عطاء رقم T47-2021.22 أجهزة مختبرات



عطاء أجهزة مختبرات

وثائق العطاء:

أ- الجزء الأول:

- (1) دعوة العطاء
- (2) الشروط والتعليمات التنظيمية للعطاء
 - (3) طريقة الدفع

ب- الجزء الثاني:

(1) جدول الكميات والمواصفات الفنية



الجزء الأول (1)

إعلان طرح عطاء رقم 747-2021.22

أجهزة مختبرات

تدعو الجامعة العربية الأمريكية الشركات المختصة الى المشاركة في العطاء المذكور أعلاه. يمكن الاستفسار أو الحصول على وثائق العطاء من دائرة اللوازم والمشتريات في الجامعة/ مبنى الدوائر الإدارية الطابق الثاني، هاتف - 2418888 -04 وثائق العطاء من دائرة اللوازم و 2510972 ونيpnp@aaup.edu مقابل مبلغ غير مسترد مقداره 100 دولار أمريكي تدفع في إحدى البنوك المعتمدة وذلك اعتباراً من يوم (الاربعاء) الموافق .15/6/2022

ملاحظات .

- 1. تقديم عرضين: فني ومالي، وسيتم دراسة العروض فنياً ومالياً لاختيار العرض المناسب.
- 2. آخر موعد لتسليم العطاءات هو في تمام الساعة الثانية من يوم(الاربعاء) 29/6/2022 ولنفس المكان.
- 3. يجب تقديم كفالة دخول عطاء 5% من قيمة العطاء على شكل كفالة بنكية أو شيك بنكي مصدق لصالح الجامعة العربية الأمريكية.
- 4. الأسعار (دولار) وتشمل جميع الضرائب بما فيها ضريبة القيمة المضافة وعلى المورد تقديم الفواتير الضريبية وشهادة خصم المصدر.
 - 5. الجامعة غير ملزمة بأقل الأسعار وبدون إبداء الأسباب.
 - 6. رسوم الاعلان على من يرسو عليه العطاء.
- 7. بإمكانكم الاطلاع على النظام الداخلي لدائرة اللوازم والمشتريات من خلال زيارة صفحة الجامعة العربية الامريكية على الانترنت. www.aaup.edu



الشروط والتعليمات التنظيمية للعطاء

(2)

- 1. على جميع المشاركين في العطاء الالتزام التام بهذه الشروط والتعليمات، وهي تعتبر جزءاً لا يتجزأ من أي أمر
 شراء أو عقد يبرم مع المشارك الفائز ما لم ينص صراحة على خلاف ذلك في أمر الشراء أو العقد.
 - 2. في هذه الشروط والتعليمات يرمز إلى "الجامعة العربية الامريكية بالاختصار (AAUP).
 - 3. يجب أن تكون الشركة المتقدمة للعطاء مسجلة رسمياً ومشتغلاً مرخصاً.
 - 4. تقدم الأسعار (دولار) شاملاً لجميع الضرائب بما في ذلك ضريبة القيمة المضافة (VAT).
- 5. يلتزم المشارك الفائز بتقديم شهادات خصم المصدر والفواتير الضريبية اللازمة وأية مستندات قانونية أخرى تغطى عملية الشراء.
- 6. يجب أن تشتمل الأسعار على جميع المصاريف المطلوبة من النقل والتركيب والتشغيل والفحص والصيائة والتدريب في المواقع المحددة في جدول المواصفات والكميات المرفق.
 - 7. يجب أن تكون الأسعار المقدمة سارية المفعول لمدة لا تقل عن (90) يوماً من تاريخ تقديم العرض.
- 8. على المشارك الفائز تقديم كفالة حسن تنفيذ خلال أسبوع من تاريخ الاتفاقية بحيث تعادل (10%) من قيمة الاتفاقية على شكل كفالة بنكية صادرة عن إحدى البنوك العاملة في فلسطين أو شيك مصدق صادر لصالح "الجامعة العربية الامريكية".
- 9.إذا تخلف المناقص الفائز عن تقديم كفالة حسن التنفيذ عن الموعد المحدد في البند السابق فإنه يحق لـ (AAUP) الغاء الاحالة.
- 10. إذا تخلف المناقص الفائز عن التوقيع على عقد التنفيذ و تسليم الكفالات والتأمينات المطلوبه منه خلال أسبوع من تاريخ قرار الاحالة، يعتبر مستنكفا عن تنفيذ العطاء ويصادر مبلغ الكفالة أو التأمين دخول العطاء بالاضافة الى ذلك يتحمل فرق السعر و/أو اي أضرار أخرى قد تلحق بالجامعة نتيجة استنكافه ويحرم من لمشاركة في عطاءات الجامعة لمدة عام.
- 11.اذا تخلف المناقص الفائز عن تنفيذ العطاء الذي احيل عليه او خالف شرطا من شروط العقد يحق للجامعة مصادرة كفالة دخول العطاء أو حسن التنفيذ أو جزء منها وتنفيذ العطاء مباشرة من الجامعة أو اية جهة تراها مناسبة بالاسعار والشروط والطريقة المناسبة ويتحمل المناقص أي فروقات بالاسعار مضاف اليها 15% من اجمالي قيمة العطاء.

12 يتحمل المناقطيل المتخلف وفي تعويض بدل اي عطل او ضرر قد يلحق بالجامعة نتيجة لذلك.

- 13. تعاد كفالة حسن التنفيذ بعد استكمال التوريد وجميع شروط العقد أو أوامر الشراء وبموجب الوثائق الأصولية اللازمة للاستلام.
 - 14. على المشاركين في العطاء ارفاق كتالوجات عن المنتج.
- 15. يلتزم من يرسو عليه العطاء بدفع غرامة تأخير بواقع (0.1%) عن كل يوم تأخير من قيمة الأعمال المنجزة عن الوقت المحدد في الاتفاقية، ويتم احتساب هذه الغرامات من الدفعات المستحقة له أو من كفالة حسن التنفيذ.
- 16. يحق لـ (AAUP) إلغاء العطاء دون إبداء الأسباب كما أن (AAUP) غير ملزمة بإحالة العطاء على أقل العروض سعراً دون إبداء الأسباب. ولها أن ترفض كل أو بعض العروض المقدمة لها دون أن يكون لأي من المشاركين الحق في الرجوع إليها بأي خسارة أو ضرر ناجم عن تقديم عرضه ولا يترتب على (AAUP) أي التزامات مادية أو غير مادية مقابل ذلك، كما يحق ل (AAUP) تجزئة العطاء بما تراه مناسبا ودون ابداء أسباب.
- 17. يلتزم من يرسو عليه العطاء بتقديم كفالة بنكية (صيانة) بقيمة (5%) من قيمة الأعمال المنجزة صالحة لمدة عام من تاريخ تسليم الأعمال.
- 18. على المشارك في العطاء تقديم عرضه على أساس المواصفات الفنية االمبينة في وثائق العطاء وبموجب الكميات المرفق.
- 19. لا يجوز للمشارك في العطاء أن يتنازل لأي طرف آخر عن كل أو جزء من أمر الشراء دون الحصول على إذن خطى من (AAUP) مع الاحتفاظ بكامل حقوق (AAUP) وفقاً لشروط أمر الشراء.
- 20. عند دراسة العروض يؤخذ بعين الاعتبار كفاءة المناقص من الناحيتين المالية والفنية وقدرته على الوفاء بالتزامات العطاء وخبرته في تقديم اللوازم المطلوبة والسمعة التجارية والتسهيلات التي يقدمها ويجوز استبعاد عرضه لنقص كل أو بعض هذه المتطلبات.
 - 21. لا تقبل العروض أو التعديلات التي ترد بعد التاريخ والموعد المحدد كآخر موعد لتقديم العروض.
 - 22. يجب تعبئة جداول المواصفات المرفقة و لن ينظر بأي عرض لا يلتزم بتعبئة الجداول.
 - ويسمح بتقديم عرضين اثنين فقط كحد اقصى لكل بند.
 - يجب تقديم عرضي الاسعار الفني والمالي بنسختين: الأولى ورقية، والأخرى الكترونية (محوسبة).
 - تقديم العرضين المالي والفني الورقيين بالظرف المختوم، مع ضرورة وضع ختم الشركة والتوقيع على كل الصفحات (للعرض المالي بالذات).



(3)

طريقة الدفيع

خلال (90) يوماً من التوريد والقبول والاستلام النهائي، مقابل تقديم الكفالات المطلوبة.



الجزء الثانى

1. جدول الكميات والمواصفات الفنية أجهزة مختبرات

No.	Product	Qty	Unit Price USD	Total Price USD
1.	Thermocycler PCR	1		The second secon
2.	Horizontal Gel Electrophoresis	1		
3.	Electrophoresis generators- Electrophoresis generator PS 30	1		
4.	Mini wide horizontal tank natural colour MT-108	1		
5.	Gel Documentation System	1		
6.	Combs adapted to the SUB10 Mini-Plus Submarine Gel Electroph	1		
7.	Combs for SUB10 mini plus submarine tank 10 x 11.5 cm	1		
8.	Electroporation Transfection	1		
9.	CUY21(NepaGene)electroporator	1		
10.	Pipette puller	1		
11.	Fiber Optic Light Source	1		
12.	veterinary anesthesia machine	1		
13.	Mini centrifuge	1		
14.	Blotting system	1		
15.	Ultra low -25 C Deep Freezer	1		
16.	fluorescent microscope	1		
17.	Hotplate and stirrer	1		
18.	Dry bath (block heater)	1		
19.	Co2 incubator	1		
20.	Inverted microscope	1		
21.	LAMINAR Air FLOW- vertical	1		
22.	Digital PCR SYSTEM (Droplet)	1		
23.	Refrigerated Bench top Centrifuge	1		
24.	Bench top Centrifuge	1		
25.	THERMAL mixer	1		
26.	FUME HOOD	1		
27.	RNA/DNA Extraction Machine	1		= 1
28.	(3 digits) balance	1	and a second	131
29.	(5 digits) analytical balance	- Annual Control	1750	1:1

30.	Oven	1	_
31.	Water bath	1	
32.	Moisture Analyzer	1	
33.	Disintegration tester	1	
34.	Hardness tester	1	
35.	Friability tester	1	
36.	Tablet coating machine	1	
37.	Eliza reader	2	
38.	desktop tablet sugar coating machine	1	
39.	digital caliper	2	
40.	Pipette (10) ml	2	
41.	Micropipette (0.2-2 µI)	6	
42.	Micropipette (2- 20 µl)	6	
43.	Micropipette (20-200µl)	6	•
44.	Micropipette (100-1000µI)	6	
	Total		

في حالة وجود استفسار يرجى تزويدنا بها من خلال البريد الالكتروني للرد عليها pnp@aaup.edu



Thermocycler PCR

- Gradient Thermo Cycler
- Temperature
 - ° range : 4 − 99 °C .
 - ° Accuracy at 90 C°: ± 0.3 °C.
 - ° Temp uniformity at 90 C° : \pm 0.4 °C .
 - ° Average Ramp rate : 2.5 °C/sec
 - ° Gradient range: 30 99 °C.
 - ° Thermal Gradient differential Range: 1 20 °C.
- · Sample block:
 - block 96 samples .
 - Sample Capacity: 0.2 ml tubes or plate or strip / 0.5 ml tubes Plate or strip .
- Display :
 - ° LCD.
 - ° Textual and graphic display.
- Program :
 - Up to 99 cycles.
 - Program Memory: more than 300.
 - With Self test program .
- Heated lid
- Interface :
 - USB
- Power: 230 V.
- With temperature calibration certificate (NIST)

* Kindly your offer will not be taken into our consideration if you don't fill below table.

Thermocycler PCR	Yes/NO	Comment
Manufacture	Manufacture	
Model	Model	
Origin	Origin	
Delivery time	weeks	
3 years Full warranty		
Gradient Thermo Cycler	Yes/No	
Temperature	Yes/No	
° range : 4 – 99 °C .	Yes/No	
° Accuracy at 90 C°: ± 0.3 °C.	Yes/No	
° Temp uniformity at 90 C° : ± 0.4 °C .	Yes/No	
° Average Ramp rate : 2.5 °C/sec .	Yes/No	
° Gradient range: 30 - 99 °C.	Yes/No.	
° Thermal Gradient differential Range : 1 - 20 °C	Yserso	111
Sample block:	Yellow	1 1

PPD-P03-R05 2 25/07/2018

° block (96) samples .	Yes/No
 Sample Capacity: 0.2 ml tubes or plate or strip / 0.5 ml tubes Plate or strip . 	Yes/No
Display:	Yes/No
° LCD.	Yes/No
° Textual and graphic display .	Yes/No
Program :	Yes/No
Up to 99 cycles .	Yes/No
Program Memory : more than 300 .	Yes/No
With Self test program .	Yes/No
Heated lid	Yes/No
Interface :	Yes/No
° USB	Yes/No
Power : 230 V	Yes/No
Complete and Ready to use with all accessories	Yes/No
User's and service manual	Yes/No
Price include installation and training in site	Yes/No
Complete and Ready to use with all accessories	Yes/No
User's and service manual .	Yes/No
Date of manufacture not less than 2022	Yes/No
10 years availability for spare parts	Yes/No



Items 2 3 4

Horizontal Gel Electrophoresis

Two Gels size: (mid size) With safety lid + (mini size)

Buffer volume: 300 or 400 ml 2 sets of combs with different size

UV transparent tray. casting trays and combs

Positive/ Negative electrode (with levelling bubble).

Power supply:

Power Supply suitable for above horizontal gel electrophoresis and basic applications

Power supply: 10 - 300 V, 4 - 500 mA, 80 W, with 1 mA steps.

with a 3 digit LED display.

Programmable methods.

Timer control.

Pause and resume function.

Safety: no load detection, sudden load change detection, overload/short circuit detection and over voltage protection.

Complete system and ready to use.

Horizontal Gel Electrophoresis with power supply	Yes/No Comment
Manufacture for power supply and tank	Manufactur
	e
Model for power supply and tank	Model
Origin: USA, Europe , japan for power supply and tank	Origin
Delivery time	Delivery
	time
3 years Full warranty	Yes/No
Complete set including :	Yes/No
Two Gels size: (mid size) With safety lid + (mini	Yes/No
size)	
Buffer volume: 300 or 400 ml	Yes/No
2x 16 samples of combs & -WELLs	
UV transparent tray.	
casting trays and combs	
Positive/ Negative electrode (with levelling bubble).	Yes/No
Power supply:	Yes/No
Power supply: 10 - 300 V, 4 - 500 mA, 80 W, with 1mA	Yes/No
steps.	
Power Supply suitable for above horizontal gel	Yes/No
electrophoresis	
with a 3 digit LED display .	Yes/No
timer control .	Yes/No
Programmable methods .	Yes/No
pause and resume function .	Yes/No
Safety: no load detection, sudden load change detection,	Yes/No
overload/short circuit detection and over voltage protection	
10 years availability for spare parts	Yes/No
Complete and Ready to use with all accessories	Yes/No
User's and service manual	Yes/No
Price include installation and training in site	Yes/No
Date of manufacture not less than 2022	Yes/N

Item 7 6 5

Gel Documentation System

Origin: USA, ERUO.

Application: research lab (molecular genetics)

Complete system : include dark room , camera filters ,Pc hardware with software (operation and

analysis)

Application: Chemiluminescence /Gel Documentation

Filters: UV transilluminator, epi white. Preset control (filter ,lens ,light ..)

High resolution camera: 4 mega pixel resolution or more

Motorized Zoom lens with broad range filter

UV and white light source: with foldable white light

Fire wire / USB connectivity with fast data transfer.

3 position filter slider.

dark room with sliding as well as open door facility

Computer: Complete system, brand name HP or Dell , as manufacture recommendation or above specification , include operation and analysis software .

Gel Documentation System	Yes/No	Comment
Manufacture	Manufacture	
Model	Model	
Origin	Origin	
Delivery time	time	
3years Full warranty	Yes/No	
Complete system : include dark room , camera , filters ,Pc	Yes/No	
hardware with software (operation and analysis)		
Application: Chemiluminescence /Gel Documentation	Yes/No	
Filters : UV transilluminator, epi white	Yes/No	
Preset control (filter ,lens ,light)	Yes/No	_
High resolution camera : 4 mega pixel resolution or more	Yes/No	le .
Motorized Zoom lens with broad range filter	Yes/No	
UV and white light source: with foldable white light	Yes/No	
Fire wire / USB connectivity with fast data transfer.	Yes/No	
3 position filter slider .	Yes/No	
dark room with sliding as well as open door facility	Yes/No	
Computer:	Yes/No	
Complete system, brand name HP or Dell , as manufacture recommendation or above specification , include operation and analysis software .	Yes/No	
10 years availability for spare parts	Yes/No	
Complete and Ready to use with all accessories	Yes/No	
User's and service manual	Yes/No	
Date of manufacture not less than 2022	Yes/No	
Price include installation and training in site	Yes/No	



Items 8 9

Electroporation Transfection

DC Pulse Waveform: Square

Impedance Measurement: Max. 35Kohm. Voltage: 0.1-99.9V (0.1V resolution)

Pulse Length: 05-99.9ms (0.01ms resolution). **pulse Interval**: 0.1-999.9ms (0.1ms resolution.

Number of Pulses: 1-99.

Safety Limit (High Current): 1.60A

Measured Output (Automatic): Voltage: 0.1-99.9V / Current: 1-999mA (1mA

resolution), .00-1.60A (0.01A resolution)

Memory: 60 programs

With suitable accessories for Mouse Embyo (In Utero)

With platinum electrodes

1-1-1 Electroporation Transfection	Yes/No	Comment
nufacture	Manufacture	
del	Model	
gin	Origin	
ivery time	Delivery time	
Pulse Waveform : Square	Yes/No '	
pedance Measurement : Max. 35Kohm .	Yes/No	
tage: 0.1-99.9V (0.1V resolution)		
se Length: 05-99.9ms (0.01ms resolution).	Yes/No	
se Interval: 0.1-999.9ms (0.1ms resolution.	Yes/Na	
mber of Pulses: 1-99.	Yes/No	
ety Limit (High Current) : 1.60A	Yes/No	
asured Output (Automatic): Voltage: 0.1-99.9V / Current	Yes/No	
999mA (1mA resolution), .00-1.60A (0.01A resolution)		
mory : 60 programs		
:h suitable accessories for Mouse Embyo (In Utero)	Yes/No	
:h platinum electrodes	Yes/No	
ears warranty	Yes/No	
years availability for spare parts	Yes/No	



Pipette puller

- · Simple and easy to use.
- Digital display
- · clearly shows status and setting.
- . Heating value can be adjusted as needs, even after the needle is placed in the heater coil.
- Acrylic cover closes securely.
- Elapsed time (from the "START" to end.) shows on display.
- Thin-walled Glass Capillary .
- With Glass Capillary with Filament (ø1mm × 90mm), Spare Heater (Kanthal wire).

Pipette puller	Yes/NO	Comment
Manufacture	Harata da re	
Model	17894-1	
<u>Origin</u>	C. gr.	
Delivery time		
3 years Full warranty		
Simple and easy to use.	1 - 1 - 1	
Digital display	j. va <u>di</u> ja	
clearly shows status and setting.	- 2	
Heating value can be adjusted as needs, even after the needle is placed in the heater coil.		
Acrylic cover closes securely.		
Elapsed time (from the "START" to end.) shows on display.	. * <u>1. m</u> g	
Thin-walled Glass Capillary .	121-c	
Price include installation training in site	Yeahra	
Date of manufacture not less than 2022	Y FEA TO	
10 years availability for spare parts	Me Walts	



Fiber Optic Light Source

- · Led light .
- Light flux 600 lm.
- Max active light guide dia :9 mm
- Wide range power supply .
- · Long lamp lime.
- Protection class II.
- Color temperature (k): 5500.

Fiber Optic Light Source	Yes/NO	Comment	
Manufacture	and Year 1, 1015		
Model	. 4		
Origin	4_ f, 1) ·		
<u>Delivery time</u>	62,55		
3 years Full warranty			
Led light .			
Light flux 600 lm.	* = = 1		
Max active light guide dia :9 mm			
Wide range power supply .	La Consta		
Long lamp lime .	Yesh		
Protection class II .	2.00		
Color temperature (k): 5500.	i sipina		
Price include installation training in site	Kero		
Date of manufacture not less than 2022	7 Su . 6		
10 years availability for spare parts	1 2 2 2 2 2 2 2		



veterinary anesthesia machine

application: veterinary Respiration mode:

Manual mode, open closed mode, semi-open & semi-closed mode

Gas source:O2,quick inlet port O2 FLUSH: 20L -75L /min adjustable Flow meter: O2,

Airway pressure gauge:-20 to 100 cmH20

APL valve : adjustable easily, Anaesthetic (isoflurane) vaporizer

Circle rebreathing system with CO2 absorption easy to view Inspiration and expiration valve

Breathing circuit pressure gauge.

With breathing bags.

With Anesthesia vaporizer.

Soda lime chamber.

Bains connector: open breath system, like Bains circuit and Jackson circuit

Portable stand with 5 castors

* Kindly your offer will not be taken into our consideration if you don't fill below table.

veterinary anesthesia machine	Yes/No	Comment
Manufacture	Manufac	
	ture	
Model	Model	
Origin:	Origin	
Delivery time	Delivery	
	time	
3 years warranty		
application : veterinary	Yes/No	
Respiration mode :	Yes/No	
Manual mode, open closed mode, semi-open & semi-closed		
mode		
Gas source:O2,quick inlet port	Yes/No	
O2 FLUSH: 20L -75L /min .	Yes/No	
adjustable Flow meter: O2.	Yes/No	
Airway pressure gauge:-20 to 100 cmH20.	Yes/No	
APL valve : adjustable easily,	Yes/No	
easy to view Inspiration and expiration valve		
With breathing bags .	Yes/No	
Anesthetic (isoflurane) vaporizer		
Circle rebreathing system with CO2 absorption		
Breathing circuit pressure gauge .	Yes/No	
With Anesthesia vaporizer .	Yes/No	
Soda lime chamber .	1,007110	
Bains connector: open breath system, like Bains circuit and	Yes/No	
Jackson circuit		
Portable stand with 5 castors	Yes/No	
Price include installation and training on site	Yes/No	
10 years availability for spare parts	Yes/No	
Price include training in site .	Yes/Maria	

PPD-P03-R05 2 25/07/2018

SIL STATE TO THE PLANT THE

Mini centrifuge

Brand name

• Origin: USA, Europe.

- Bench top mini centrifuge with rotor.
- · Brushless motor.
- Microprocessor control.
- Max Speed: 2900 RPM.
- Rotor: for microtube and PCR strip with adaptors 0.2, 0.4, 0.5 ml

With Control and display.

Mini Centrifuge	Yes/No	Comment
Manufacture	Manufacture	
Model	Model	
Origin	Origin	
Delivery time	weeks	
3 years Full warranty 10 years availability for spare parts	Yes/No	
For research lab genetics lab .	Yes/No	
Bench top mini centrifuge with rotor.	Yes/No	
Brushless motor.	Yes/No	
Microprocessor control .	Yes/No	
Max Speed: 2900 RPM.	Yes/No	
Rotor: for microtube and PCR strip with adaptors 0.2, 0.4	Yes/No	
,0.5 ml		
With Control and display .	Yes/No	
Date of manufacture the machine not less than 2022		
User's and service manual .	Yes/No	
Price include installation and team training in site.	Yes/No	



Blotting system

brand name

origin: USA.EURO
application: Genetics lab
turbo blotting machine
fast and semi dry.
4 or 2 gels / run.

Mini protein cell include:

tank, lid power cables, electrode, mini cell buffer dam, 2 casting stands with 4 casting frames, 10 glass plates for

1mm thick gels, 10 well combs. gel release.

Suitable higher power supply for above

Gel holder cassettes

Foam pad . Cooling unit .

Blotting system	Yes/No	Comment
Manufacture	Manufactur	
	е	
Model	Model	
Origin	Origin	
Delivery time	Delivery	
	time	
3 years warranty		
Country of origin (USA , Europe)	Yes/No	
application : Genetics lab	Yes/No	
turbo blotting machine	Yes/No	
fast and semi dry .	Yes/No	
4 or 2 gels / run .	Yes/No	
Mini protein cell include :	Yes/No	
tank, lid power cables, electrode, mini cell buffer	Yes/No	
dam , 2 casting stands with 4 casting frames , 10		
glass plates for		
1mm thick gels , 10 well combs . gel release .	Yes/No	
Suitable higher power supply for above	Yes/No	
Gel holder cassettes	Yes/No	
10 years availability for spare parts	Yes/No	
Price include training in site .	Yes/No	



<u>Ultra low -25 C Deep Freezer</u>

Brand name.

Origin: USA, Europe, japan. research lab application.

Up right .

Volume: minimum 280 L. Temp . Range : up to – 25 C .

Shelves: 4.

Fully programmable microprocessor controlled with membrane keypad and eye level control panel.

Refrigeration system: Cascade.

Minimum temperature recovery time after door opening.

lockable outer doors and lids are required .

Insulation: polyurethane 90 mm thick.

Must have battery back-up.

Audible and visible alarms for: Temp. High/low, filter check, self-diagnostics, door check, power failure, battery low.

Smart diagnostics for simpler system diagnostics

Must use CFC-FREE, HCFC-FREE non-flammable refrigerants.

Freezers should have heated air vent and front panel air filter.

Heavy-duty lockable castors.

Noise Level: less than 50 (db).

Freezer must have standard quality test requirements and Electrical safety certified.

With Temperature calibration certificate.

Warranty: 3 years

Optional:

Analog chart recorder, data logger or Digital chart recorder, 7 different size Racks with Boxes. Online UPS enough for 30 minutes working time.

* Kindly your offer will not be taken into our consideration if you don't fill below table.

Ultra low -25 C Deep Freezer	Yes/No	Comment
Manufacture	Manufacture	
Model	Model	
Origin: USA, Europe , japan	Origin	
Delivery time	Delivery time	
Warranty: 3 years	Yes/No	
For biology and research lab application.	Yes/No	
Up right .	Yes/No	
Volume: minimum 280 L.	Yes/No	
Temp . Range : up to - 25 C .	Yes/No	
Shelves: 4.	Yos/No	
Fully programmable microprocessor controlled with membrane keypad and eye level control panel.	Yes/No	and the second s
: Temp monitoring system	1	

PPD-P03-R05 2 25/07/2018

COUTEM OF AND FURCHASING

Refrigeration system : Cascade .	Yes/No	
Minimum temperature recovery time after door opening.	YesiNo	
lockable outer doors and lids are required.	Yes/No	
Insulation: polyurethane 135 mm thick.	Yes/No	_
Must have battery back-up .	Yes/No	
Audible and visible alarms for: Temp. High/low, filter check, self-diagnostics,	Yes/No	
door check, power failure, battery low,		
Smart diagnostics for simpler system diagnostics	Yes/No ·	
Must use CFC-FREE, HCFC-FREE non-flammable refrigerants.	Yes/No	
Freezers should have heated air vent and front panel air filter.	Yes/No	
Heavy-duty lockable castors.	Yes/No	
Noise Level: less than 50 (db).	Yes/No	
Freezer must have standard quality test requirements and Electrical safety certified.	Yes/No	
With Temperature calibration certificate.	Yes/No	
Optional:	Yes/Ne	
Analog chart recorder , data logger or Digital chart recorder , 7		
different size Racks with Boxes.Online UPS enough for 30		
minutes working time		
10 years availability for spare parts	Yes/No	
Complete and Ready to use with all accessories	YosiNo	
User's and service manual	Yes/No	
Price include installation training in site	Yes/No	



fluorescent microscope

- Brand name
- Origin: USA. EUOPE. JAPAN.
 Example for our requested (axio 1mager 2 z2) الموديل المعروض لدلالة فقط على المطلوب
- Warranty: 3 years.
- Application : Research lab (for Cell structure visualization, and neurobiology research
)
- · For continuous operation .
- System include a device used to take optical sections in fluorescence imaging. These individual section images are combined together in a computer to create a 3D image.
- · Optics :
- Techniques : bright field , phase contrast, DIC , florescence
- · Magnification 6 objectives :
- (10x, 20x, 63x) water, (40x, ,63 x 100x) oil, (20x) air.
- Wide field focusing eyepieces of view with soft rubber eyecups and should be suitable for spectacle wearers. The eyepieces should be suitable for graticule insertion.
- The filters should be easily insert able and removable in the reflector. Filters: Fluorescence for UV, Blue and Green (49 dpi , 38he(GFP) , 43HE(DSRED) , 50 (CY5) , 46HE(YFP) , 47 (CFP) , 62HE(BFP/GFP) AND DIC TRANS LIGHT)
- Dust cover, all wires, cords, connectors and standard accessories needed for proper functioning of the microscope
- Optional :
 - With camera
 - With Computer and analysis software (Image capture and analysis software that shall be compatible with the camera)

fluorescent microscope	Yes/No	Comments
Manufacture	Manufacture	
Model	Model	
Origin	Origin	
Delivery time	weeks	
Warranty 3 years	Yes/No	
Application : Research lab (for Cell structure visualization, and neurobiology research)	Yes/No	
For continuous operation .	Yes/No	
 System include a device used to take optical sections in fluorescence imaging. These individual section images are combined together in a computer to create a 3D image. 	Yes/No	

• Optics :	Yes/No
Techniques: bright field, phase contrast, DIC, florescence	Yes/No
Magnification 6 objectives :	Yes/No
• (10x , 20x , 63x) water , (40x , , 63 x 100x) oil , (20x) air .	Yee/No
Wide field focusing eyepieces of view with soft rubber eyecups and should be suitable for spectacle wearers. The eyepieces should be suitable for graticule insertion.	
The filters should be easily insert able and removable in the reflector. Filters: Fluorescence for UV, Blue and Green (49 dpi, 38he(GFP), 43HE(DSRED), 50 (CY5), 46HE(YFP), 47 (CFP), 62HE(BFP/GFP) AND DIC TRANS LIGHT)	
Dust cover, all wires, cords, connectors and standard accessories needed for proper functioning of the microscope	
User's and service manual .	Yes/No
Price include installation and team training in site.	Yes/No
References in Palestine are Attached to your offer	Yes/No
10 years availability for spare parts	Yes/No

Hotplate and stirrer

- Brand name, USA, Europe origin
- Ceramic heating plate with protective layer and resistance to chemicals.
- Speed Control: up to 1000 rpm.
- Stirring volume: up to 2L.
- With stirrer bar
- Temperature range control: up to 200 C°.
- Optional price:

Heater, thermostat, stirring bar and stir bar removal, holding rod.

Hot plate stirrer	Yes/No	Comment
<u>Manufacture</u>	Manufactur	
7	9	
<u>Model</u>	Model	
<u>Origin</u>	Origin	
Delivery time	Delivery	
	time	
3 years Full warranty		
Ceramic heating plate with protective layer and	Yes/No	
resistance to chemicals.		
Speed Control: up to 1000 rpm.	Yes/No	
With stirrer bar		
Stirring volume: up to 2 L.	Yes/No	
Temperature range control : up to 200 C°.	Yes/No	
Power : 220 V .	Yes/No	
Price include installation training in site	Yes/No	
User's and service manual	Yes/No	
10 years availability for spare parts	Yes/No	
Optional price :		
Heater , Thermostat , stirring bar and stir bar	Price	
<u>removal</u> , holding rod .		



Dry bath (block heater)

Brand name .

Origin: USA, Europe. For Heating test tubes. Temp range: 5-165 C.

Uniform and stable temp: at 37c (0.1C±).

Fast recovery time for pre-setting heating temp .

Heater 300 W.

Heating aluminium block size : (24x0.5 ml, 24x1.5 ml, $24 \times 2 \text{ ml}$ + $1 \times 96 \text{ well}$).

Dry bath (block heater)	Yes/No	Comment
<u>Manufacture</u>	Manufacture	
<u>Model</u>	Model	
<u>Origin</u>	Origin	
Delivery time	time	
3 years Full warranty	Yes/No	
For Heating test tubes and well plates	Yes/No	90
Temp range : 5 – 165 C .	Yes/No	
Uniform and stable temp: at 37c (0.1C±).	Yes/No	
Fast recovery time for pre-setting heating temp .	Yes/No	
Heater 300 W .	Yes/No	
Heating aluminium block size : (24x0.5 ml ,	Yes/No	
24x1.5 ml , 24 x 2 ml + 1x 96 well) .		
10 years availability for spare parts	Yes/No	
Complete and Ready to use with all accessories	Yes/No	



Co2 incubator

Brand name

Origin: USA, Europe .

- Temperature range up to +50 °C.
- Volume: 85 L.
- microprocessor temperature controller.
- Electronic microprocessor temperature controller with auto-diagnostic system, programme stored in case of power failure.
- Direct heat: 4 sides with back and door heating
- Clear display: temperature, time, CO2, humidity and set-up values.
- Fast recovery time for co2 and Temperature .
- Heater: 1500 W
- Alarm function: visual and acoustic alarm, integral fault diagnostics (CO2 and empty gas cylinder, humidity, temperature, open door)
- · Serialization program.
- Co2 measurement : Infrared
- Internal chamber: made of stainless steel
- · Inner glass door .
- 3 perforated stainless steel shelves
- Over temperature protection.
- digital electronic CO2-control with automatic zeroing (infrared measuring system) adjustment range 0 10 % resolution 0.1 %
- opening hole for test purposes
- · interruption CO2-supply if the door opening.
- With one CO2 gas cylinder and hose and meter gauge and co2 filter.
- · With one point calibration certificate.

Co2 incubator	Yes/No	Comment
Manufacture	Man.	
Model	Model	
Origin	Origin	
Delivery time	D.time	
Country of origin (USA , Europe)	Yes/No	
3 years warranty	Yes/No	
Temperature range up to +50 °C .	Yes/No	
Volume: 85 L.	Yes/No	
microprocessor temperature controller.	Yes/No	
Electronic microprocessor temperature controller with auto- diagnostic system , programme stored in case of power failure .	Yes/No	,
Direct heat: 4 sides with back and door heating	Yes/No	
Clear display: temperature, time, CO2, humidity and set-up values.	Yes/No	
Fast recovery time for co2 and Temperature .	Yes/No	
Heater : 1500 W	Yes/No	

Alama for the second and a second and the second an	A # 10 - 2 16 11 -	
Alarm function: visual and acoustic alarm, integral fault	Yes/No	Į.
diagnostics (CO2 and empty gas cylinder , humidity ,		
temperature, open door)		
Serialization program.	Yes/No	
Co2 measurement : Infrared	Yes/No	
Internal chamber: made of stainless steel	Yes/No-	
Inner glass door .	Yes/No	
3 perforated stainless steel shelves	Yes/No	
Over temperature protection .	Yes/No	
digital electronic CO2-control with automatic zeroing (infrared	Yes/No	
measuring system) - adjustment range 0 - 10 % - resolution 0.1		
%		
opening hole for test purposes	Yes/No	
interruption CO2-supply if the door opening.	Yes/No	
With one CO2 gas cylinder and hose and meter gauge and	Yes/No	
co2 filter .		
With one point calibration certificate.	Yes/No	
Price include installation and team training in site.	Yes/No	
User's and service manual .	Yes/No	
10 years availability for spare parts	Yes/No	



Inverted microscope

- Brand name :
- Origin: USA. EUOPE. JAPAN.
- Warranty: 3 years.
- Techniques : bright field(BF) , phase contrast (PC).
- · inverted microscope for application cell culture
- Nosepiece : Quintuple (5 position), backward-facing type.
- Stand :

Rugged and sturdy inverted stand with low positioned co-axial coarse and fine focus knobs with harmonic focusing. Transparent stage .

• Eyepiece :

2 eyepiece wide filed (with pointer) of view with swivelling eyepiece, rubber eyecups and should be suitable for spectacle wearers . tube is inclined 35° and rotate 360°

- Optics : infinity optic
- Objective: Magnification: 4x, 10x, 20x, 40x,...
- Light: LED Illumination.
- Stage: large stage plate with mechanical stage of traveling. Should be able to hold petridish and specimen slides or cell culture flask.
- Dust cover, all wires, cords, connectors and standard accessories needed for proper functioning of the microscope.

Inverted microscope	Yes/No	Comments
Manufacture	Manufacture	
Model	Model	
Origin	Origin	
Delivery time	weeks	
Warranty 3 years	Yeshlo	
Techniques : bright field(BF) , phase contrast (PC) ,	Yes/No	
Nosepiece: Quintuple (5 position), backward-facing type. Stand: Rugged and sturdy inverted stand with low positioned co-axial coarse and fine focus knobs with harmonic focusing. Transparent stage.	Yes/No	
Eyepiece: 2 eyepiece wide filed (with pointer) of view with swivelling eyepiece, rubber eyecups and should be suitable for spectacle wearers. tube is inclined 35° and rotate 360°.	Yes/No	
Optics : infinity optic	Yes/No	
Objective: Magnification (BF & PC): 4x, 10x, 20x, 40x.	Yes/No	
Light: LED Illumination .	Yes/No	
Stage: large stage plate with mechanical stage of traveling. Should be able to hold petridish and specimen slides or cell culture flask	Yes/No	
Dust cover, all wires, cords, connectors and standard accessories needed for proper functioning of the microscope	Yes/No	,
User's and service manual .	Yes/No	
Price include installation and team training in site.	Yes/No	
References in Palestine are Attached to your offer	Yes/No	71
10 years availability for spare parts	A CONTRACTOR OF THE PARTY OF TH	18

LAMINAR Air FLOW- vertical

Brand name .

Horizontal Laminar flow .

Bench top: 4 ft.

light: UV, Florescence

Air purification : HEPA Filter: 99.999 %efficiency removal 0.3 um practical I removal .

Microprocessor control.

Filter replacement warning signal .

Door should be made of tempered safety glass sliding door or wind screen .

Noise Level: <40 dBA.

Working area: good quality Stainless steel.

Blower: High efficient centrifugal type with lifetime lubricated bearings.

Electric supply: 220 V.

Should have FDA CE approval.

Optional price:

HEPA filter . UV Lamp . Fluorescent light .

LAMINAR Air FLOW- vertical .	Yes/No	Comment
Manufacture	Manufactur	
	e	
Model	Model	
Delivery time	Delivery	
	time	
Warranty : 3 years.	Yes/No	
For research lab application.	Yės/No	
Brand name .	Yes/No	
Origin:	Yes/No	,
Horizontal Laminar flow .	Yes/No	
Bench top: 4 ft.	Yes/No	
light: UV, Florescence	Yes/No _	
Air purification: HEPA Filter: 99.999 %efficiency removal	Yes/No	
0.3 um practical I removal .		
Microprocessor control .	Yes/No	
Filter replacement warning signal .	Yes/No	
Door should be made of tempered safety glass sliding door	Yes/No	
or wind screen .		
Noise Level : <40 dBA.	Yes/No	
Working area: good quality Stainless steel.	Yes/No	
Blower : High efficient centrifugal type with lifetime	Yes/No	
lubricated bearings		
Electric supply: 220 V.	Yes/No	
Should have FDA CE approval .	Yes/No	
Optional price :	Yes/No	
HEPA filter . UV Lamp Fluorescent light .	Yes/No	
Complete and Ready to use with all accessories	Yes/No	
User's and service manual	Yes/No	
Price include installation training in site	Yes/No	



Digital PCR SYSTEM (Droplet)

Application : Cancer biomarker studies , Pathogen detection , Next generation sequencing ,

Gene expression analysis.

Complete system include: Droplet generator, Droplet reader, computer, software.

SAMPLE SIZE : 20 UL .

Generator Capacity: 1 to 8 samples / cartridge

Droplet per 20 UI: 20,000

reader capacity: up 96 samples.

Sample illumination: LED

Detection channel: FAM (EvaGreen), HEX (VIC)

Precision: 10%.

Easy access to previously stored sealing protocols.

* Kindly your offer will not be taken into our consideration if you don't fill below table.

Droplet digital PCR SYSTEM	Yes/NO	Comment
Manufacture		
Model	1 = 1 = 2 = 22 =	
Origin	* Mid of	
Delivery time	1.45 - 7.	
3 years Full warranty		
Application: Cancer biomarker studies, Pathogen detection, Next		
generation sequencing , Gene expression analysis.		
Complete system include: Droplet generator, Droplet reader,		
computer, software.		
SAMPLE SIZE : 20 UL	J. 14	
Generator Capacity: 1 to 8 samples / cartridge		
Droplet per 20 UI: 20,000	5 1 ·	
reader capacity: up 96 samples.	d * 1	
Sample illumination : LED		
Detection channel : FAM (EvaGreen), HEX (VIC)	V	
Precision: 10%		
User's and service manual with ALL software password.		
Price include installation and full training in site	100 T	
Date of manufacture not less than 2022	1111111	
10 years availability for spare parts	14, 14	
Complete with all standers accessories to working with our application		



PPD-P03-R05 2 25/07/2018

Page 29 of 52

Refrigerated Bench top Centrifuge

- Origin: USA, Europe.
- Bench top centrifuge.
- · Brushless motor .
- Microprocessor control.
- MAX Speed: RFC 24,500 & 16000 RPM
- Low noise 50 Db
- Temp Range: (-10 or lower to +35 °C)
- Control and display:
 - Digital display.
 - actual values are indicated During centrifugation .
 - Speed display: RPM or RFC, input in steps of 10
 - RCF key: indication during run.
 - Time (t/min): 1 99 min, continuous operation.
 - Impulse key: For short centrifugation operations .
- Centrifuge chamber: stainless steel.
- Rotor:
 - Easily exchangeable rotors: Capability to use two types of rotor swing out and fixed.
 - swing out rotor with buckets and lids .
 - samples: tubes (15 ml: 17 Ø x20 l ml, speed 4500 RPM), angle 90.
 falcon tube ((50 ml: 29 Ø x115 l ml, speed 4500 RPM), angle 90.
 EDTA tube ((4-7 ml: 13 Ø x100 l ml, speed 5000 RPM), angle 90.

EDTA tube ((1,5 -5 ml : 13 Ø x75l ml , speed 5000 RPM) , angle 90 . EDTA tube ((8,5 -10 ml : 16 Ø x100l ml , speed 4000 RPM) , angle 90 . PCR PALTE 96 + PCR strip 0.2 (speed 4000 RPM) , angle 90 .

- <u>fixed rotor</u> with bio-containment <u>samples</u>: tubes (2 ml , 1.5 ml ,0.5 ml , 0.2 ml , speed 14000 RPM) , angle 45 . tubes (, 15 ml ,50 ml, speed 9000 RPM) , angle 35 .

Safety:

- Emergency lid lock release.
- Imbalance switch-off.
- Automatic rotor recognition .
- Power: 220 V.

* Kindly your offer will not be taken into our consideration if you don't fill below table.

Refrigerated Centrifuge	Yes/No	Comment
<u>Manufacture</u>	Manufacture	
<u>Model</u>	Model	
<u>Origin</u>	Origin	
Delivery time	weeks	
3 years warranty 10 years availability for spare parts	Yes/No	
Bench top centrifuge with rotor.	Yes/No	
Brushless motor .	Yes/No	
Microprocessor control.	Yes/No	
• Temp Range: (-10 or lower to +35 °C)	and a	
Low noise 50 Db		

PPD-P03-R05 2 25/07/2018

Page **30** of **52**

Tender 47/2021.22 AAUP

Control and diaplay	VaniNa	
Control and display: Digital display.	Yes/No	
. ,	Yes/No	
- actual values are indicated During centrifugation .		
- Speed display: RPM or RFC, input in steps of 10		
- RCF key: indication during run.		
- Time (t/min): 1 - 99 min, continuous operation.		
- Impulse key: For short centrifugation operations		*
Centrifuge chamber : stainless steel .	Yes/No ·	
Rotor:	Yes/No	
Easily exchangeable rotors: Capability to use two types	Yes/No	
of rotor swing out or fixed angle		
• Rotor :	Yes/No	
- Easily exchangeable rotors : Capability to use two		
types of <u>rotor swing out and fixed</u> .		
 swing out rotor with buckets and lids. 		
- samples: tubes (15 ml:17 Ø x20 l ml, speed		
4500 RPM), angle 90 .		
falcon tube((50 ml:29 Ø x115 T		
ml,speed 4500 RPM),angle 90.		
EDTA tube ((4-7 ml : 13 Ø		
x100 l ml ,speed 5000 RPM),angle 90.		
EDTA tube ((1,5-5 ml : 13 Ø x75l		
ml,speed 5000 RPM),angle 90.		
EDTA tube ((8,5 -10 ml : 16 Ø		
x100l ml ,speed 4000 RPM), angle 90 .		
PCR PALTE 96 + PCR strip 0.2 (
speed 4000 RPM),angle 90.		
- <u>fixed rotor</u> with bio-containment		
samples: tubes (2 ml , 1.5 ml ,0.5 ml , 0.2 ml ,		
speed 14000 RPM), angle 45.		
tubes (, 15 ml , 50 ml, speed 9000 RPM		
) , angle 35 .		
Safety:	Yes/No	
Emergency lid lock release, Imbalance switch-off.	Yes/No	
Automatic rotor recognition	Yes/No	
Power: 220 V.	Yes/No	
User's and service manual .	Yes/No	
Price include installation and team training in site.	Yes/No	
	Yes/No	



Bench top Centrifuge

- Origin: USA, Europe.
- · Bench top centrifuge.
- Brushless motor .
- Microprocessor control.
- MAX Speed: RFC 24,500 & 16000 RPM
- Low noise 50 Db
- Control and display:
 - Digital display.
 - actual values are indicated During centrifugation .
 - Speed display: RPM or RFC, input in steps of 10
 - RCF key: indication during run.
 - Time (t/min): 1 99 min, continuous operation.
 - Impulse key: For short centrifugation operations .
- Centrifuge chamber: stainless steel.
- Rotor:
 - Easily exchangeable rotors: Capability to use two types of rotor swing out and fixed.
 - swing out rotor with buckets and lids .
 - samples: tubes (15 ml: 17 Ø x20 l ml, speed 4500 RPM), angle 90. falcon tube ((50 ml: 29 Ø x115 l ml, speed 4500 RPM), angle 90. EDTA tube ((4-7 ml: 13 Ø x100 l ml, speed 5000 RPM), angle 90.

EDTA tube ((1,5-5 ml : 13 Ø x75l ml , speed 5000 RPM) , angle 90 . EDTA tube ((8,5-10 ml : 16 Ø x100l ml , speed 4000 RPM) , angle 90 . PCR PALTE 96 + PCR strip 0.2 (speed 4000 RPM) , angle 90 .

- <u>fixed rotor</u> with bio-containment

<u>samples</u>: tubes (2 ml , 1.5 ml ,0.5 ml , 0.2 ml , speed 14000 RPM) , angle 45 . tubes (, 15 ml , 50 ml, speed 9000 RPM) , angle 35 .

- Safety:
 - Emergency lid lock release.
 - Imbalance switch-off.
 - Automatic rotor recognition .
- Power: 220 V.

* Kindly your offer will not be taken into our consideration if you don't fill below table.

<u>Centrifuge</u>	Yes/No	Comment
Manufacture	Manufacture	
Model	Model	
Origin	Origin	
Delivery time	weeks	
3 years warranty 10 years availability for spare parts	Yes/No	
Bench top centrifuge with rotor.	Yes/No	
Brushless motor .	Yes/No	1
Microprocessor control .	Yes/New war	
Low noise 50 Db		

PPD-P03-R05 2 25/07/2018

of 52

Control and display:	Yes/No
Tomas and anopias	The state of the s
 Digital display . actual values are indicated During centrifugation . Speed display : RPM or RFC , input in steps of 10 RCF key : indication during run . Time (t/min) : 1 - 99 min, continuous operation . Impulse key: For short centrifugation operations 	<u>Yes/No</u>
Centrifuge chamber : stainless steel .	Yes/No
Rotor:	Yes/No
Easily exchangeable rotors : Capability to use two types	Yes/No
of rotor swing out or fixed angle	
 Rotor: Easily exchangeable rotors: Capability to use two types of rotor swing out and fixed. swing out rotor with buckets and lids. samples: tubes (15 ml: 17 Ø x20 l ml, speed 4500 RPM), angle 90.	Yes/No .
Safety:	Yes/No
Emergency lid lock release, Imbalance switch-off.	Yes/No
Automatic rotor recognition	Yes/No
Power: 220 V.	Yes/No
User's and service manual .	Yes/No
Price include installation and team training in site.	Yes/No
	Yes/No



THERMAL mixer

- Brand name,
- Speed control: 200 2000 RPM.
- With heating up to 90 c.
- up to 18 to 30 samples at one time
- Universal rack for small and large samples

Thermal Mixer	Yes/NO	Comment
Manufacture	Manufacture	
Model	Model	
Origin	Origin	
Delivery time	Delivery	
	time	
3 years warranty 10 years availability for spare parts	Yes/No	
Application: General Laboratory.	Yes/No	
Vortex Tube .	Yes/No	
Speed control : 200 – 2000 RPM	Yes/No	
With heating up to 90c.		
up to 18 to 30 samples at one time		
up to 18 samples at one time		
Universal rack	Yes/No	
User's and service manual .	Yes/No	
Price include installation and team training in site.	Yes/No	



FUME HOOD

Hood to protects operator from toxic vapor and environment

Application: for general purpose

Working Dimension: 130 cm w x 70 cm D x 75 cm.

Microprocessor control (clear and simple).

Construction and Material: Made of good quality melamine board with good acid and alkali resistance function anti chemical, heat, fire resistant and easily cleanable.

Sash: Tempered Glass vertical rising sash.

electric socket (water proof).

Lighting: fluorescent lamp with vapor proof design.

Blower:

Built-in centrifugal blower for air suction in hood, silent (compatible with above specification).

Adjustable air speed.

casters and locking for casters.

with stand.

Noise level: 55 db.

Ducting:

- · chemical resistant rigid and flexible ductwork from hood to exhaust .
- التركيب يشمل توصيل الجهاز لتمديدات وتركيب في الموقع المخصص

Fume hood	Yes/No	Comment
Manufacture	Manufacture	
Model	Model	
Origin	Origin	
Delivery time	Delivery	
5years Full warranty	time Yes/No	
	Yes/No	
Hood to protects operator from toxic vapor and environment	Yes/No	
Application: for general purpose chemistry Lab		
Working Dimension: 130 cm w x 70 cm D x 75 cm.	Yes/No	
Microprocessor control (clear and simple).	Yes/No	
Construction and Material: Made of good quality melamine board with	Yes/No	
good acid and alkali resistance function anti chemical, heat, fire		
resistant and easily cleanable.		
Sash: Tempered Glass vertical rising sash.	Yes/No	
electric socket (water proof).	Yes/No	
Lighting: fluorescent lamp with vapor proof design.	Yes/No	
Blower:	Yes/No	
Built-in centrifugal blower for air suction in hood, silent (compatible with		
above specification).		
Adjustable air speed.		
casters and locking for casters.	Yes/No	
with stand.	Yes/No	
Noise level: 55 db.	The state of the s	
Ducting:	18	
chemical resistant rigid and flexible ductwork from hood to exhaust.	Yes/No	

Fume hood	Yes/No	Comment
Price include installation and training in site	Yes/No	
Complete and Ready to use with all accessories	Yes/No	
User's and service manual	Yes/No	
Date of manufacture not less than 2022	Yes/No	



RNA/DNA Extraction Machine

- Brand name :
- Origin: USA. EUOPE. JAPAN.
- · Warranty: 3 years.
- General: automated isolation of genomic DNA/RNA from whole human blood and other tissue samples
- Application: Research lab (for purification of nucleic acid DNA & RNA, proteins and cell isolation) (DNA and RNA extraction. High throughput sequencing compatible quality and)
- volumes: 200 ul up to 3mL
- Yield: Samples need to be in a least concentration of 50 ng/Microliter, total DNA yield of (at least) 3micrograms of dsDNA.
- Extraction of nucleic acids from 3-4 cc blood.
- Sample can be whole blood, fetal amniontic fluid cells. Fetal cells circulating in maternal blood.
- Automated magnetic bead technology.
- · fast with minimum hands-on time .
- with UV lamp: easy and effective decontamination (exposure time up to 16 hour)
- Open System for reagents and consumables and software for extraction protocols
- Heating and cooling function
- 6 or 12 samples /run or Number of Samples: Up to 16
- Procedure duration: up to 60 minutes.
- Magnetic Particle Size: ca. 1µm
- Can run two purification methods simultaneously
- Internal memory for ca. 150 protocols
- · Complete system with all reagents and consumables and software .

Kindly your offer will not be taken into our consideration if you don't fill below table

RNA/DNA Extraction Machine	Yes/No	Comments
Manufacture	Manufacture	
Model	Model	
Origin	Örigin	
Delivery time	weeks	
Warranty 3 years	Yes/No	
 General: automated isolation of genomic DNA/RNA from whole human blood and other tissue samples 	Yes/No	
 Application: Research lab (for purification of nucleic acid DNA & RNA, proteins and cell isolation) (DNA and RNA extraction. High throughput sequencing compatible quality and) 	Yas/No	
 volumes : 200 ul up to 3mL 	Yes/No	
 Yield: Samples need to be in a least concentration of 50 ng/Microliter, total DNA yield of (at least) 3micrograms of dsDNA. 	Yes/No	
 Extraction of nucleic acids from 3-4 cc blood . 	Yes/No	
 Sample can be whole blood, fetal amniontic fluid cells. Fetal cells circulating in maternal blood. 	Value of the same	

Automated magnetic bead technology .	Yes/No	1 1
fast with minimum hands-on time .		
with UV lamp: easy and effective decontamination (exposure time up to 16 hour)		
Open System for reagents and consumables and software for extraction protocols		
Heating and cooling function		
6 or 12 samples /run or Number of Samples: Up to 16		
Procedure duration: up to 60 minutes.	Yes/No	
Magnetic Particle Size: ca. 1µm		
Can run two purification methods simultaneously		
Internal memory for ca. 150 protocols		
Complete system with all reagents and consumables and software.		
User's and service manual .	Yes/No	
Complete system with all reagents and consumables and software		
Price include installation and team training in site.	Yes/No	
References in Palestine are Attached to your offer	Yes/No	
10 years availability for spare parts	Yes/No	

(3 digits) balance

- Brand name, USA, Europe origin
- Capacity: 2000 g.
- with cover
- Readability: 0.001 g
- 3 digits.
- Internal & external calibration.
- Display :LCD .
- Unit: Milligram, Gram,
- · Stainless steel pan (resistance to chemicals).
- Power : AC adaptor (supplied) or internal rechargeable battery.
- · With calibration certificate.
 - * Kindly your offer will not be taken into our consideration if you don't fill below table.

3 digit Balance	Yes/No	Comment
Manufacture	Manufactur	
Model	Model	
Origin	Origin	
Delivery time	Delivery time	
3 years warranty and 10 years availability for spare parts	Yes/No	
Capacity: 2000 g.	Yes/No	
Readability :0.001g	Yes/No	
With cover	Yes/No	
Display :LCD.	Yes/No	
Unit: Milligram, Gram,	Yes/No	
Stainless steel pan (resistance to chemicals).	Yes/No	
Internal & external calibration	Yes/No	
Power : AC adaptor (supplied) or internal rechargeable battery	Yes/No	
With calibration certificate .	Yes/No	



(5 digits) analytical balance

- Brand name, USA, Europe origin
- Capacity: 200 g.
- with cover
- Readability: 0. 01 mg
- Display :LCD .
- Unit: Milligram, Gram,
- Stainless steel pan (resistance to chemicals).
- Power: AC adaptor (supplied) or internal rechargeable battery.
- With calibration certificate.
 - * Kindly your offer will not be taken into our consideration if you don't fill below table.

Balance	Yes/No	Comment
Manufacture	Manufactur	,
	0	
Model	Model	
Origin	Origin	
Delivery time	Delivery	
	time	
Capacity: 200 g.	Yes/No	
Readability :0.01mg	YesiNo	
With cover		
Display :LCD.	Yes/No	
Unit: Milligram, Gram,	Yes/No	
Stainless steel pan (resistance to chemicals).	Yes/No	
Power : AC adaptor (supplied) or internal	Yes/No	
rechargeable battery		
Power: 220 V.	Yes/No Yes/No	
Internal and external calibration .		
User's and service manual .	Yes/No	
Price include installation and team training in site.	Yes/No	
3 years warranty	Yes/No	
10 years availability for spare parts	Yes/No	



Oven

- Country of origin (USA, Europe) and brand name.
- For general lab application.
- Volume 6 0L .
- **Temp range:** up to 30 250°C
- Microprocessor control .
- Adjustable parameter : time , temp ,
- Self diagnostics for any error.
- · Clear display.
- resolution of display for setpoint and actual temperature values 0.1°C.
- · program stored in case of power failure .
- Fast recovery time .
- · Two shelves .
- Access port.
- With temperature calibration certificate.

* Kindly your offer will not be taken into our consideration if you don't fill below table.

<u>Oven</u>	Yes/No	Comment
Manufacture	Manufactu	
	1,6	
Model	Model	
Origin	Origin	
Delivery time	Delivery	
	time	
3years Full warranty	Yes/No	
Country of origin (USA , Europe)	Yes/No	
For general lab application .	Yes/No	
Microprocessor control	Yes/No	
Volume 60 L .	Yes/No	
Temp range: up to 30 - 250°C	Yes/No	
Adjustable parameter : time , temp ,	Yes/No	
Self diagnostics for any error .	Yes/No	
<u>Clear display</u> .	Yes/No	
resolution of display for set point and actual temperature values 0.1°C	Yes/No	
<u> -</u>		
program stored in case of power failure .	Yes/No	
Fast recovery time .	Yes/No	п
Two shelves .	Yes/No	
With temperature calibration certificate .	Yes/No	
10 years availability for spare parts	Yes/No	
Service and user manual	Yes/No	
Price should be include installation and team training in site	Yes/No	



Waterbath

- Volume
 - 22 L.
- Control
 - microprocessor temperature controller.
 - temperature range : 5 99.9 C□ with resolution 0.1.
 - integrated autodiagnostic system with fault indicator .
 - Digital display (LED) of set and actual temperature
 - Digital timer: 1 min- 95,59 hours.
 - Low water Sensor.
- Visual alarm :
 - Over temperature protection .
- Interior chamber (I x w x h):
 - (570 x 439 x 293 mm) easy-to-clean interior material 1.4301 (ASTM 304), corrosion resistant and proof
- Exterior chamber (I x w x h) :
 - (349 x 288 x 219 mm) made of stainless steel, material no. 1.4301 (ASTM 304), laser-welded ,corrosion resistance and proof.
- Electrical :
 - 220V 50Hz.
 - 1700 W .
- With separate thermometer.
- · Cleaning spray for cleaning water bath .

Optional price:

- Shaking device :
 - Shaker control: shaking frequency (10 130 oscillation /min) with holder and clips for flasks volume (100, 200, 300,500 ml).
 - Adaptor for flacon tube 15, 50 ml.
 - 1.5 ml, 2ml.
 - Stainless steel Gable cover
 - Kindly your offer will not be taken into our consideration if you don't fill below table.

WATER BATH	Yes/No	Comments
Manufacture		
Model		
Origin		
3years Full warranty		
Control :		
temperature range : $5 - 99.9$ C with resolution 0.1.		
integrated autodiagnostic system with fault indicator .		
Digital display (LED) of set and actual temperature		The state of the s
Digital timer: 1 min- 95,59 hours.		and the same of th

WATER BATH	Yes/No	Comments
Low water Sensor .		
Visual alarm :		
Over temperature protection .		
Interior chamber (I x w x h):		
(570 x 439 x 293 mm) easy-to-clean interior material 1.4301 (ASTM 304), corrosion resistant and proof		
Exterior chamber (I x w x h) :		
(349 x 288 x 219 mm) made of stainless steel, material no. 1.4301 (ASTM 304), laser-welded ,corrosion resistance and proof. VOLUME: 22 LITERS		
Electrical : 220V 50Hz. 1700 W		
Cleaning spray for cleaning water bath .		
Optional price :		-
Shaking device :		
 Shaker control: shaking frequency (10 – 130 oscillation /min) with holder and clips for flasks volume (100, 200, 300,500 ml). Adaptor for flacon tube 15, 50 ml. 1.5 ml, 2ml. 		
- Stainless steel Gable cover		
Price include installation and training in site		



Moisture Analyzer

Brand name, USA, Europe, japan

Capacity: 50 g.

instrument to determine relative humidity for small samples of different materials.

Readability: 0.001 g

Heating: halogen

Chamber size: height 8mm and 90 diameter.

Temp: up to 150 C.

Products: up to 900 products

Weighings: up to 1000 weighings and Programs: up to 100 drying programs

results obtained using the moisture analyzer: % M, % D, % R, g (sample mass after drying.)

Display :clear display

With calibration certificate.

Kindly your offer will not be taken into our consideration if you don't fill below table.

Moisture Analyzer	Yes/No	Comment
Manufacture	Manufacture	NEW WORLD SERVICE STREET, STRE
Model	Model	
Origin	Origin	
Delivery time	Delivery time	
3 years warranty and 10 years availability for spare parts	Yes/No	
Capacity : 50 g .	Yes/No	
instrument to determine relative humidity for small samples of different materials .	Yes/No	
Readability :0. 001 g	Yes/No	
Heating : halogen	Yes/No	
Chamber size : height 8mm and 90 diameter .	Yes/No	
Temp: up to 150 C.	Yes/Hone	
Products : up to 900 products	res/Nd	

Moisture Analyzer	Yes/No	Comment
Weighings : up to 1000 weighings and Programs : up to 100 drying programs	Yes/No	
results obtained using the moisture analyzer : % M , % D , % R , g (sample mass after drying.)	Yes/No	
Display :clear display	Yes/No	
System must be ready and complete with all accessories to work completely	Yes/No	
Price include instillation and training in site with service and user manual		



Disintegration tester

Brand name, USA, Europe,, japan

tablet disintegration tester.

With two basket clear view and easy clean.

With magnetic coupling.

Quick lock

Digital heart (independent).

Low liquid level.

Basket for large tablet and capsules

Clear display

Complete system with all accessories and ready to use

• Kindly your offer will not be taken into our consideration if you don't fill below table.

Disintegration tester	Yes/No	Commen
Manufacture	Manufacture	
Model	Model	
Origin	Origin	
Delivery time	Delivery time	
3 years warranty and 10 years availability for spare parts	Yes/No	
tablet disintegration tester.	Yes/No	
With two basket clear view and easy clean.	Yes/No	
With magnetic coupling.	Yes/No	
Quick lock	Yes/No	
Digital heart (independent).	Yes/No	
Low liquid level .	Yes/No	
Basket for large tablet and capsules	Yes/No	
Clear display	Yes/No	
System must be ready and complete with all accessories to work completely	Yes/No	
Price include instillation and training in site with service and user manual		

Hardness tester

Brand name, USA, Europe, japan.

• Measuring units: breaking force, diameter, weight, thickness

• LCD display.

• Test ranges: 0 - 300 N

• TEST SPEED: 1-40mm/min 0.5±.

• Units: N, KP,Kg, lbs.

• Batch test.

Kindly your offer will not be taken into our consideration if you don't fill below table.

Hardens tester	Yes/No	Comment
Manufacture	Manufacture	
Model	Model	
Origin	Origin	
Delivery time	Delivery time	
3 years warranty and 10 years availability for spare parts	Yes/No	7
LCD display.	Yes/No	
Test ranges: 0 – 300 N	Yes/No	
TEST SPEED: 1-40mm/min 0.5±.	Ves/No	
Units: N, KP,Kg, lbs.	Yes/No	
Batch test.	Yes/Nu	
	Yes/No	
System must be ready and complete with all accessories to work completely	Yes/No	
Price include instillation and training in site with service and user manual		



Friability tester

Brand name, USA, Europe, japan.

Stainless steel housing .

Two drums variable speed: 20 55 rpm.

With abrasion drum

Kindly your offer will not be taken into our consideration if you don't fill below table.

Friability tester	Yes/No	Comment
Manufacture	Manufacture	
Model	Model	
Origin	Origin.	
Delivery time	Delivery	
3 years warranty and 10 years availability for spare parts	Yes/No	
Stainless steel housing.	Yes/No	
Two drums variable speed .	Yes/No	
With abrasion drum	Yes/No	
System must be ready and complete with all accessories to work completely	Yes/No	
Price include instillation and training in site with service and user manual		



Tablet coating machine

Brand name, USA, Europe, japan.

Complete system: stainless steel coating pan, drive unit control panel, gun, hot air blower, stainless steel body.

Simple

Pan size: from 8" to 18"

Operation temp: $30 - 80 \, \text{C}$.

Spray System: The Unit with pump with Spray gun & silicon tubing's. Pump is a Variable speed pump provided with digital controller for precise speed regulation and rapid priming facility

Speed Controller: Speed Knob between 0 to 40 RPM.

Exhaust Blower: 0.25 HP, 1400 RPM Single Phase Exhaust blower of 100 CFM for venting the heat out of your working room.

Gear Motor Drive 0.5 HP, 1370 RPM/ Sturdy & Maintenance Free Geared Motor for continuous operation.

Single phase

• Kindly your offer will not be taken into our consideration if you don't fill below table.

	Yes/No	Comment
Manufacture	Manufacture	
Model	Model	
Origin	Origin	
Delivery time	Delivery time	
3 years warranty and 10 years availability for spare parts	Yes/No	
Complete system: stainless steel coating pan, drive unit control panel, gun, hot air blower, stainless steel body.	Yes/No	
Simple	Yes/No	A LANGE OF THE PARTY OF THE PAR
Pan size: from 8" to 18 "	Yes/News	
Speed: 0-30 rpm.	Yes/No (

	Yes/No	Comment
Single phase	Yes/No	
	Yes/No	
	Yes/No .	
	Yes/No	
Price include instillation and training in site with service and user manual		



Eliza reader

Origin: USA, Europe, Japan. Wide applicability for life science.

Open system Kits.

Could Read: 6-, 12-, 24-, 48-, 96 plate "flat, U, V" shape.

Fully Automatic ELIZA microplate reader

Read Method: Endpoint, kinetic,.

5 filters or more: between 400 - 750 nm

Light Sources: Xenon Flash Lamp . Photometric range : 0.0–3.5 OD .

Accuracy : \leq 1.0% or 0.010 from 0.000–3.000 OD at 490 nm . Linearity : \leq 1.0% from 0.0–2.0 OD ; \leq 2.0% from 0.0–3.0 OD .

Resolution: 0.001 OD.

LAMP : Tungsten halogen lamp With wide and clear Display

Built in Shaker.

Incubator with programmable temperature control.

Fast microplate reading speeds Stores 100 user entered tests.

Complete and Ready to use with all accessories.

Optional Price:

Lamp.

• Kindly your offer will not be taken into our consideration if you don't fill below table.

Eliza reader	Yes/No	Comment
Manufacture	Manufacture	
Model	Model	
Origin: USA , Europe .		
Delivery time	weeks	
3years Full warranty	Yes/No	
Brand name	Yes/No	
Open system Kits .	Yes/No	
Could Read: 6-, 12-, 24-, 48-, 96 plate '' flat , U , V " shape .	Yes/No	
Fully Automatic ELIZA microplate reader	Yes/No	
Read Method: Endpoint, kinetic,.	Yes/No	
5 filters or more : between 400 - 750 nm	Yes/No	
Light Sources: Xenon Flash Lamp .	Yes/No	
Photometric range : 0.0–3.5 OD .	Yes/No	
Accuracy : ≤1.0% or 0.010 from 0.000–3.000 OD at 490 nm .	Yes/No	
Linearity : ≤1.0% from 0.0–2.0 OD ; ≤2.0% from 0.0–3.0 OD .	Your Market	12
Resolution: 0.001 OD.	Teg Too	1.5

LAMP : Tungsten halogen lamp	Yes/No	
With wide and clear Display	Yes/No	
Built in Shaker .	Yes/No	
Incubator with programmable temperature control.	Yes/No	
Fast microplate reading speeds	Yes/No	
Stores 100 user entered tests.	Yes/No	
Complete and Ready to use with all accessories.	Yes/No	
Price include installation and training in site .	Yes/No	
User's and service manual .	Yes/No	
Optional Price:	Yes/No	
Lamp.	Yes/No	
	Yes/No	

