

دائرة اللوازم والمشتريات

عطاء رقم (2019.20/ T20) توفير أجهزة ولوازم للعلاج الطبيعي

2019-2020



عطاء توفير أجهزة ولوازم للعلاج الطبيعي

وثائق العطاء:

- أ- الجزء الأول:
- (1) دعوة العطاء
- (2) الشروط والتعليمات التنظيمية للعطاء
 - (3) طريقة الدفع

ب- الجزء الثاني:

(1) جدول الكميات والمواصفات الفنية



الجزء الأول (1)

إعلان طرح عطاء رقم (T20/2019.20)

عطاء توفير أجهزة ولوازم للعلاج الطبيعي

تدعو الجامعة العربية الأمريكية الشركات المختصة الى المشاركة في العطاء المذكور أعلاه. يمكن الاستفسار أو الحصول على وثائق العطاء من دائرة اللوازم والمشتريات في الجامعة/ مبنى الدوائر الإدارية الطابق الثاني، هاتف- 2418888 -04 وثائق العطاء من دائرة اللوازم و2510972 وثي الكترونيpnp@aaup.edu مقابل مبلغ غير مسترد مقداره 50 دو لار أمريكي تدفع في إحدى البنوك المعتمدة وذلك اعتباراً من يوم (الاثنين) الموافق 2019/11/11

ملاحظات:

- 1. تقديم عرضين: فني ومالي، وسيتم دراسة العروض فنيا ومالياً لاختيار العرض المناسب.
- 2. آخر موعد لتسليم العطاءات هو في تمام الساعة الثانية من يوم (الاثنين) 2019/11/25 ولنفس المكان.
- يجب تقديم كفالة دخول عطاء 5% من قيمة العطاء على شكل كفالة بنكية أو شيك بنكي مصدق لصالح الجامعة العربية الأمريكية.
- 4. الأسعار بالشيكل وتشمل جميع الضرائب بما فيها ضريبة القيمة المضافة وعلى المورد تقديم الفواتير الضريبية وشهادة خصم المصدر.
 - 5. الجامعة غير ملزمة بأقل الأسعار وبدون إبداء الأسباب.
 - 6. رسوم الاعلان على من يرسو عليه العطاء.



(2)

الشروط والتعليمات التنظيمية للعطاء

- على جميع المشاركين في العطاء الالتزام التام بهذه الشروط والتعليمات، وهي تعتبر جزءاً لا يتجزأ من أي أمر شراء أو عقد يبرم مع المشارك الفائز ما لم ينص صراحة على خلاف ذلك في أمر الشراء أو العقد.
 - 2. في هذه الشروط والتعليمات يرمز إلى "الجامعة العربية الامريكية بالاختصار (AAUP).
 - 3. يجب أن تكون الشركة المتقدمة للعطاء مسجلة رسمياً ومشتغلاً مرخصاً.
 - 4. تقدم الأسعار (بالشيكل) شاملاً لجميع الضرائب بما في ذلك ضريبة القيمة المضافة (VAT).
- 5. يلتزم المشارك الفائز بتقديم شهادات خصم المصدر والفواتير الضريبية اللازمة وأية مستندات قانونية أخرى تغطى عملية الشراء.
- والتدريب في المواقع المحددة في جدول المواصفات والكميات المرفق.
 - 7. يجب أن تكون الأسعار المقدمة سارية المفعول لمدة لا تقل عن (90) يوماً من تاريخ تقديم العرض.
- 8. على المشارك الفائز تقديم كفالة حسن تنفيذ خلال أسبوع من تاريخ الاتفاقية بحيث تعادل (10%) من قيمة الاتفاقية على شكل كفالة بنكية صادرة عن إحدى البنوك العاملة في فلسطين أو شيك مصدق صادر لصالح "الجامعة العربية الامريكية".
- 9. إذا تخلف المناقص الفائز عن تقديم كفالة حسن التنفيذ عن الموعد المحدد في البند السابق فإنه يحق لـ (AAUP) الغاء الإحالة.
- 10. تعاد كفالة حسن التنفيذ بعد استكمال التوريد وجميع شروط العقد أو أوامر الشراء وبموجب الوثائق الأصولية اللازمة للاستلام.
 - 11. على المشاركين في العطاء ارفاق كتالوجات عن المنتج.
- 12. يلتزم من يرسو عليه العطاء بدفع غرامة تأخير بواقع (0.1%) عن كل يوم تأخير من قيمة الأعمال المنجزة عن الوقت المحدد في الاتفاقية، ويتم احتساب هذه الغرامات من الدفعات المستحقة له أو من كفالة حسن التنفيذ.
- 13. يحق لـ (AAUP) إلغاء العطاء دون إبداء الأسباب كما أن (AAUP) غير ملزمة بإحالة العطاء على أقل العروض سعراً دون إبداء الأسباب. ولها أن ترفض كل أو بعض العروض المقدمة لها دون أن يكون لأي من المشاركين الحق في الرجوع إليها بأي خسارة أو ضرر ناجم عن تقديم عرضه ولا يترتب على (AAUP) أي التزامات مادية أو غير مادية مقابل ذلك.



- 14. يلتزم من يرسو عليه العطاء بتقديم كفالة بنكية (صيانة) بقيمة (5%) من قيمة الأعمال المنجزة صالحة لمدة عام من تاريخ تسليم الأعمال.
- 15. على المشارك في العطاء تقديم عرضه على أساس المواصفات الفنية االمبينة في وثائق العطاء وبموجب الكميات المحددة في جدول الكميات المرفق.
- 16. لا يجوز للمشارك في العطاء أن يتنازل لأي طرف آخر عن كل أو جزء من أمر الشراء دون الحصول على إذن خطي من (AAUP) مع الاحتفاظ بكامل حقوق (AAUP) وفقاً لشروط أمر الشراء.
- 17. عند دراسة العروض يؤخذ بعين الاعتبار كفاءة المناقص من الناحيتين المالية والفنية وقدرته على الوفاء بالتزامات العطاء وخبرته في تقديم اللوازم المطلوبة والسمعة التجارية والتسهيلات التي يقدمها ويجوز استبعاد عرضه لنقص كل أو بعض هذه المتطلبات.
 - 18. لا تقبل العروض أو التعديلات التي ترد بعد التاريخ والموعد المحدد كآخر موعد لتقديم العروض.
 - 19. يجب تعبئة جداول المواصفات المرفقة و لن ينظر بأي عرض لا يلتزم بتعبئة الجداول.
 - ٠٠ ويسمح بتقديم عرضان اثنان فقط كحد اقصى لكل بند
 - يجب تقديم عرضي الاسعار الفني والمالي بنسختين: الأولى ورقية، والأخرى الكترونية (محوسبة).
 - تقديم العرضين المالي والفني الورقيين بالظرف المختوم، مع ضرورة وضع ختم الشركة والتوقيع على كل الصفحات (العرض المالي بالذات)



(3)

طريقة الدفيع

خلال (90) يوماً من التوريد والقبول والاستلام النهائي، مقابل تقديم الكفالات المطلوبة.



الجزء الثانى

1. جــدول الكميـات والمواصفات الفنية عطاء توفير أجهزة ولوازم للعلاج الطبيعي

الرقم	الصنف /العدد	العدد	الوحدة	السعر/شيكل	المجموع/شيكل	الوصف	صورة توضيحية
1	Paraffin Wax Manufacture:	10	عدد			شمع بارافين يستخدم في جهاز البرافين قابل للذوبان	
2	Thera-Band Manufacture: Model: Origin:	3	لفة			مطاط بالوان مختلفة ومقاومات مختلفة تستخدم في التمارين	■TheraBand ■
3	Hot Cold Discrimination Kit Manufacture: Model: Origin:	6	Each			أداة تستخدم لفحص الإحساس عند المرضى وتستخدم الحار والبارد	Everth Healthcare
4	Mobile Mirrors (4 castors) Manufacture: Model: Origin: 1007+003 + العلاج المائي	3	Each			مراية عجلات متحركة اطارها خشبي	
5	Transfer Board Manufacture: Model: Origin:	3	Each			لوحات خشبية لنقل المريض من السرير الى الكرسي والعكس	T D
6	Hydrocollator (Cold Pack) مرفق المواصفات الفنية وجدول المطابقة	1	Each			وحدة الكمادات الباردة تحتوي بداخلها عدة احجام من الكمادات الباردة بهدف العلاج	
7	Gloves Latex (large) Manufacture: Model: Origin:	5	علبة			work glovesفازات طبية	EN
8	Infrared Lamp مرفق المواصفات الفنية وجدول المطابقة	1	عدد				With Chinese-made 271W holls

الرقم	الصنف /العدد	العدد	الوحدة	السعر/شيكل	المجموع/شيكل	الوصف	صورة توضيحية
9	Electrical traction back, neck with all accessory مرفق المواصفات الفنية وجدول المطابقة غرفة 003 وهو بدل القديم المتلف	1	عدد			جهاز يستخدم لشد الفقرات العنقية والظهرية	
10	corner stair درج خشب Manufacture : Model : Origin :	1	عدد			درج خشبي يتخدم لتدريب المريض على المشي والطلوع والنزول من الدرج ويستخدم في التاهيل	
11	Electronic scale ميزان	5	عدد			ميزان بسيط لقياس الوزن فقط دون الطول	
12	مخدة اكريلان Manufacture : Model : Origin :	30	عدد			مخدة طبية بمواصفات لتفادي الم الرقبة والعمود الفقري	
13	נرج سرير Manufacture : Model : Origin :	10	عدد			درج احادي وثناني يستخدم لصعود المريض على السرير ممكن ان يكون من الخشب او الحديد	Burloy
14	Dual Inclinometer Manufacture: Model: Origin:	10	Each			instruments that can be used in daily practice to measure active cervical range of motion (ACROM) in patients with non- specific neck pain.	
15	Wash Bottles 250 ml Manufacture: Model: Origin:	20	عدد		and and an order to the	علب بلاستيكية سعة 250مل	

الرقد	الصنف /العدد	العدد	الوحدة	السعر/شيكل	المجموع/شيكل	الوصف	صورة توضيحية
16	Wooden Wall Ladder Manufacture : Model : Origin : السعر يشمل التركيب	3	عدد			سلم حانط مصنوع من الخشب الصلب باطوال مختلفة يستخدم للتمارين العلاجية للكبار والصغار	
17	Massage Tools (set) Manufacture: Model: Origin:	2	عدد			أدوات خشب مختلفة الاشكال تستخدم في المساج والتدليك	0 Y A 4 Y O 8 8
18	Mobilization Belt Manufacture: Model: Origin:	15	عدد			Mulligan Mobilisation Belt for patient stabilization High-quality nylon belt Offers one- handed adjustment and a quick release Can be used to mobilize patient during manual therapy techniques Measures 8' in length	6
19	Goniometer Stainless Steal (set) Manufacture: Model: Origin:	2	عدد			أدوات تستخدم لقياس المدى الحركي للمفاصل مصنوعة من الستانليس باشكال مختلفة	
20	Hand sanitizer for protection agaist germ 500 ml Manufacture: Model: Origin:	15	عدد			علب بلاستيكية تحتوي على مادة معقمة للايدي	© =
21	Therapeutic balls (different sizes)	2	عدد	<u> </u>		These innovative and adaptive	R

الرقم	الصنف /العدد	العدد	الوحدة	السعر/شيكل	المجموع/شيكل	الوصف	صورة توضيحية
	Manufacture : Model : Origin :					exercise balls are ideal for immediate rehabilitation following an injury, as well as meditative daily stretching & muscle strengthening. Thera-Band Exercise Balls come with two plugs, an inflation adapter and a helpful poster illustrating 24 exercises. Three different sizes of balls are offered to accommodate user height variations.	
22	pediatric walkers and k-walker Manufacture: Model: Origin:	2	336	a bet annual and a state of the		Lightweight and height adjustable, this walker can also be easily folded away when not in use. Featuring easy to grip vinyl handgrips designed just for smaller hands, this walker also has foot piece silencers for noise-free use.	

دائية اللوازم والمشتوبات

الرقم	الصنف /العدد	العدد	الوحدة	السعر/شيكل	المجموع/شيكل	الوصف	صورة توضيحية
						including forearm supports, a seating unit and hip location pads	
23	Pediatric rolator Manufacture: Model: Origin:	2	عدد			Lightweight, sturdy aluminum frame with proportionally large width to length provides wide stance and greater stability. The Children's Walker has three inch wheels to help overcome normal floor obstacles, cracks, carpet or outside play areas, and the large rear rubber tips provide effective braking capability. The height adjustable design accommodates a wide range of user heights.	
24	Balance beam kit Manufacture: Model: Origin:	2	عدد			this tool to build the balance challenge. individual beams can be used as hurdles to jump over, or stand them up to use as bowling pins to knock down. The Balance Beam Kit is great to allow each child to construct their own unique obstacle course, addressing motor planning, coordination and of course balance work	BALANCE BEAM KIT S INOSAUR PT APPROVED!




الرقم	الصنف /العدد	العدد	الوحدة	السعر /شيكل	المجموع/شيكل	الوصف	صورة توضيحية
25	Rollers Manufacture: Model: Origin:	2	عدد			Functional tool to practice transitional movements, such as sitting to stand, and sitting to quadruped. Great to practice weight bearing through hands and feet in preparation for crawling.	OAM ROLLER = DINGSAVR PT A
26	Therapeutic exercise mattress Manufacture: Model: Origin:	2	عدد			Thin and lightweight with a 1800mm length, the Airex® Fitline 180 Mat can be used outdoors indoors or even in water. Ideal for exercise or therapy, this closed cell foam mat provides with a hard-wearing, slip resistant surface	
27	Wedges cushion Manufacture: Model: Origin:	3	عدد			tool for dynamic sitting activities, promoting postural control, trunk strength and dynamic balance	WEDGE CUSHION DINGSAUR PT ABERO
28	Wall bar for child Manufacture: Model: Origin:	2	عدد			Metallic or wood child wall bar	



الرقم	الصنف /العدد	العدد	الوحدة	السعر/شيكل	المجموع/شيكل	الوصف	صورة توضيحية
29	Stepping floor bar Manufacture: Model: Origin:	1	عدد			Metallic parallel bar with hand rails and separation board with obstacle course	● Contacts
30	Balance board Manufacture: Model: Origin:	2	عدد			The balance board engages intrinsic foot musculature, promotes lower extremity strength, alignment and balance reactions	BALANCE BOARD S
31	Bobath bed Manufacture: Model: Origin:	1	عدد			Electrical height adjustment with footswitch and central locking wheels	
32	Pediatric Playing table for child Manufacture: Model: Origin:	1	عدد			This table and chair set include a square-shaped table, whiteboard, chalkboard, two small chairs, two fabric storage bins, and two fabric holders. The wooden construction ensures optimal strength and durability	
33	Pediatric chair Manufacture: Model: Origin:	4	عدد			The Chair Cube encourages social interaction as well as gross and fine motor skill development. The two seating heights offer plenty of options for kids ages 3 and up.	
34	Balance walking mat Manufacture: Model: Origin:	1	عدد	624	and a section with an art transport to an angular	This flat Figure 8 Walking Path offers an extra coordination dimension because your clients will	3

الرقم	الصنف /العدد	العدد	الوحدة	السعر/شيكل	المجموع/شيكل	الوصف	صورة توضيحية
						have to follow a curved one-footed path. The board is made of birch and is hinged in the center so it can be folded and easily stored.	
35	Various Pediatric toys Manufacture: Model: Origin:	2	عدد			Playing and stimulation toys	BALLOONS = DINOSAUR PT APPROVED! POP TOORS INOSAUR PT APPROVED! BUSINES - DINOSAUR PT APPROVED! PLAYER BALLS = DINOSAUR PT APPROVED!
36	Pediatric Wooden toys Manufacture: Model: Origin:	2	عدد			Different sizes, colors and shapes from those wooden toys	
37	Pediatric bicycles Manufacture: Model: Origin:	2	عدد			Tricycle with hand rail ,holding rail and 3 wheels as wide base for children	WIND COMPTICES
38	Straps (for upper and lower limbs) Manufacture: Model: Origin:	10	335	(A. 1)		Different straps are used for correct the abnormal positions in the upper and lower limbs	

الرقم	الصنف /العدد	العدد	الوحدة	السعر/شيكل	المجموع/شيكل	الوصف	صورة توضيحية
39	Pediatric parallel bar Manufacture: Model: Origin:	1	عدد	السعر إشيدن	المجموع إشيدن	The Parallel Bars can be used to help teach children with mobility problems how to better bear weight on the legs, stretch the lower limbs, and relax the hips, leading to better walking and balancing skills. To suit clients with different needs, the height of the pediatric by myself component can be adjusted and they can be locked onto the Walk-On Base at various widths. Product Benefits: Enhances balance Builds leg strength Assists walking skills	صوره توصیحیه
40	pediatric positioning seat Manufacture: Model: Origin:	1	عدد			Pediatric seating system pictured with the Long Leg Sitting Board that can be used permanently or intermittently for very young children or those needing an extra stretch to the thighs or calves	
41	Standing frame Manufacture: Model: Origin:	1	عدد			This encourages upright standing by working all muscles, and reducing a hypotonic attitude. The height is adjustable with table and support for the knees and trunk	



الرقم	الصنف /العدد	العدد	الوحدة	السعر/شيكل	المجموع/شيكل	الوصف	صورة توضيحية
42	Pediatric stair Manufacture: Model: Origin:	1	عدد			This Pediatric Step and Ramp combination has been specially designed to help with gait balance and cognitive training. On one side there are three, 12cm (4.75") high steps and on the other side there is a 1:8 gradient ramp, each are fitted with a non-slip surface to maximize the child safety when in use	
43	Peanut ball Manufacture: Model: Origin:	2	عدد			The peanut ball is a great tool to practice tall kneeling, quadruped position, and prone propping	PEANUT BALL S DINOSAUR PT APPROVED
44	incline ramp Manufacture : Model : Origin :	1	عدد			This is a piece of equipment to encourage tummy time, crawling, climbing, and standing balance. Helpful for practicing early motor milestones with young children, encouraging dynamic movement. The Incline Ramp is also excellent for use in treatment of toe walking and flat feet	Trecline Panp -
45	Floor scooter Manufacture: Model: Origin:	4	acc	(E)	Hamad or that from the square bands to	Helpful tool to use both seated and prone. Addressing overall strength, motor planning, coordination and body awareness. Connect a few scooters together to	PLOOP SCOOTER DINOSAUR RY APPROVED!

idaginally sailed de la

الرقم	الصنف /العدد	العدد	الوحدة	السعر/شيكل	المجموع/شيكل	الوصف	صورة توضيحية
						create new fun activities with or without peers	
46	Foot prints Manufacture: Model: Origin:	2	عدد			Offering visual and tactile cues for foot placement to incorporate in many activities, addressing alignment, motor planning, body awareness and promoting activation of the muscles needed to maintain stable base of support	TACTILE FOOTPRIME
47	Stepping stones Manufacture: Model: Origin:	10	عدد			for creating obstacle courses, practicing step up/step downs, working on balance and motor planning. For older children allow them to set up the Stepping Stones in the configuration of their choice to encourage active participation from the planning stage to the execution.	STEPPING STON
48	Dynamic disc Manufacture: Model: Origin:	2	عدد			exercise balance device, DynaDisc® is an inflatable pillow capable of withstanding 2000lbs/sq. inch. It's smooth on one side and textured with pegs on the other. Use DynaDisc® seated to give you benefits similar to a gymball, standing to increase balance and muscular resistance while doing exercises, or lying to improve	

A mark to the first the product of

الرقم	الصنف /العدد	العدد	الوحدة	السعر/شيكل	المجموع/شيكل	الوصف	صورة توضيحية
						strength and core stabilization for floor exercises	
49	Bucket bridge Manufacture: Model: Origin:	1	عدد			for navigating, addressing alignment and lower extremity strength. Can be modified to increase step length, practice stepping in and out, stepping and jumping up and down	BUCKET BRIDGE = DINOSAUR PT APPROVED
50	Shape bean bags Manufacture: Model: Origin:	2	عدد			These Shape Bean Bags are a tool for children of all ages, providing helpful sensory input, allowing child to identify colors, shapes and for our older children reading the words on each bean bag. Bean bags allow us to practice hand eye coordination, object manipulation and grading of force and movement in throwing skills	Shape Bean Bags = Dinosaur PT Approved San Bogs
51	Dead weights Manufacture: Model: Origin:	2	عدد			Supplied as a pair, these 0.5kg (1 lb) Strap on Exercise Weights are ideal physiotherapy aids for exercising the upper and lower extremities	0.5kgs
52	Building blocks Manufacture: Model: Origin:	2	عدد	The fall and the first state of the fall and	Company of the control of the contro	Plastic building blocks and colorful toys with different size	

Company of the state of the sta

الرقم	الصنف /العدد	العدد	الوحدة	السعر/شيكل	المجموع/شيكل	الوصف	صورة توضيحية
53	Music toys Manufacture: Model: Origin:	4	عدد			the kit is supplied in a strong and sturdy storage box and includes enough instruments for groups of up to twenty people encouraging inclusion in a group environment.	
54	rings set Manufacture : Model : Origin :	1	عدد			18Pcs Kid's Stacking Rings Ring Toss Cast Throw Circle Set Toys	
55	Therapeutic bands Manufacture: Model: Origin:	5	عدد			The elastic bands that are used for strength training and muscle building, as well as for physical therapy. In general, they consist of a rubber band about 10 to 15 inches long. Some come with handles at each end which can be gripped with your hands or placed at the end of your feet, although there are some models without handles	
56	Therapeutic swing Manufacture : Model : Origin : السعر يشمل التركيب	1	عدد		المنافذة الم	The freestanding therapeutic swing offers your facility all the benefits of a ceiling suspension system to accommodate a wide range of therapy modules for vestibular, neurodevelopmental	the state of the s
			1			, and sensory integration therapy	

indiquilalland

الرقم	الصنف /العدد	العدد	الوحدة	السعر/شيكل	المجموع/شيكل	الوصف	صورة توضيحية
57	Basketball set Manufacture: Model: Origin:	1	عدد			This set Includes: basketball hoop & kid sized ball	¥I:
58	Trampoline Manufacture: Model: Origin:	1	عدد			Using the trampoline with a bar to start can help children develop independent jumping skills, recruiting upper body support until they are ready to rely on lower body strength to propel vertically	AOD. THOMIG. WWW
59	Alcohol 70% Manufacture: Model: Origin:	10	لتر			كحول معقم تركيز 70%	arear near
60	Ultrasound unit مرفق المواصفات الفنية وجدول المطابقة	1	Each				
61	Laser Therapy مرفق المواصفات الفنية وجدول المطابقة	1	Each				ENDOLASER 422 S Name of a common amount of the co
62	Electrical Stimulation unit. مرفق المواصفات الفنية وجدول المطابقة	1	Each		and a more and a fire and a second	The state of the s	

Control AND PURCHASING

الرقم	الصنف /العدد	العدد	الوحدة	السعر/شيكل	المجموع/شيكل	الوصف	صورة توضيحية
	العلاج المائي						
63	Hot Pack USA Manufacture: Model: Origin:	12	عدد				
64	Vaccum Cup Electrodes, Vaccum Sponges Manufacture: Model: Origin:	4	عدد				000
65	Walking ladder Manufacture : Model : Origin : السعر يشمل التركيب	3	عدد				at their port
	المجموع الكلي						

• يشترط بجميع الأجهزة والأدوات والمستهلكات ان تكون المواد المستخدمة والمصنعة بها آمنه عند الاستخدام مع المرضى ولايوجد بها حواف او ضعف تسبب في أي ملاحظات ضرر

ويشترط التركيب والتدريب لجميع الأجهزة عند التسليم.

يجب تعبئة المرفقات (جدول المطابقة) بشكل واضح عند التسعير .



Hydrocollator

Origin: USA, EURO.

Volume: 50 L.

• Temperature range: 71-75 C.

Thermal Cut-Out Temp: 80 - 85 C.

Heat up Time (to 71°C): 8 hours/day.

· fiberglass insulation to prevent heat loss.

· high grade stainless steel.

Simple to fill with water (No plumbing required).

· With water drain for easy cleaning.

• Included in your unit price set of (3 HotPac and Foam Filled Terry).

. Optional price and cat No:

heater.

* يرجى التكرم بتعبئة الجدول بالأسفل بالإضافة الى تقديم عرض تقنى خاص بالشركة

Hydrocollator	Yes/No	Comment
Manufacture	Manufacture	
Model	Model	
Origin	Origin	
Delivery time	Delivery time	
Two years warranty and 10 years availability for spare parts	Yes/No	
Applications: rehabilitation.	Yes/No	
Volume: 50 L.	Yes/No	
Temperature range: 71–75° C.	Yes/No	
Thermal Cut-Out Temp: 80- 85° C.	Yes/No	
Heat up Time (to 71°C): 8 hours/day.	Yes/No	
fiberglass insulation to prevent heat loss.	Yes/No	
high grade stainless steel.	Yes/No	
With water drain for easy cleaning.	Yes/No	
Your Price included with U.P: set of (3 HotPac and Foam Filled Terry).	Yes/No	
Simple to fill with water (No plumbing required) .	Yes/No	
Optional (Price & cat No):		
Heater price with cat No	FRICE	
Price include user and service manual	Yes/No	
Price include installation and training in site .	Yes/No	
Date of manufacturing not less than 2019		

تلتزم الشركة بتقديم كافة ملحقات والاضافات الازمة لعمل الجهاز تسعيرها بشكل واضح وابداء اية ملاحظات ضرورية للجهاز أو المكان قبل التوريد

The state of the s

Pediatric Platform Mounted Parallel Bars

Application: physiotherapy. Origin: USA, EURO.

- · hardwood plywood base with ulta-durable
- used to help teach children with mobility problems how to be better bear weight on legs stretch lower limbs and relax the hips and better walking and balancing.
- · Parallel bars mounted stainless steel handrails for easy use and durability
- Safety tread at both ends of the platform
- Heavy duty, square, telescoping, steel uprights
- · Easy height changes with spring-loaded plungers and numbered height indicator strips
- Second height knob removes side to side play
- · End bumpers on all handrails
 - Kindly your offer will not be taken into our consideration if you don't fill below table.

Platform Mounted Parallel Bars	Yes/No	Comment
Manufacture	Manufacture	
Model	Model	
Origin	Origin	
Delivery time	Delivery	
Two years warranty and 10 years availability for spare parts	Yas/No	
hardwood plywood base with ulta-durable	Yes/No	
used to help teach children with mobility problems how to be better bear weight on legs stretch lower limbs and relax the hips and better walking and balancing		
Parallel bars mounted diameter stainless steel handrails for easy use and durability	Yes/Na	
Safety tread at both ends of the platform	Yes/No	
Heavy duty, square, telescoping, steel uprights	Yes/No	
Easy height changes with spring-loaded plungers and numbered height indicator strips		
Second height knob removes side to side play		
End bumpers on all handrails		
Price include installation and training in site.		

تلتزم الشركة بتقديم كافة ملحقات والاضافات الازمة لعمل الجهاز أو تصعيرها بشكل واضح وابداء اية ملاحظات ضرورية للجهاز أو المكان قبل التوريد

مادة العالم والمستويات

Traction machine

Origin: USA, EURO.

- · For lumbar and cervical.
- · User defined rest , treat time , hold .
- · Real time digital treatment tracking ..
- Therapy: continuous, intermittent, harmonized, progressive.
- Traction force: 0-90 kg.
- Force time: / rest time: up to 98s.
- Treat time: 0 to 90 min.
- Safety: audible signal (patient interrupt switch activated, end of treatment)
- With all accessories: cervical traction attachment, pelvic attachment, knee bolsters, Flexion Stools, traction cord with metal fitting, patient interrupt switch, cervical harness separated bar, include in main price.
- Optional price :
 Examination bed for retraction

• يرجى التكرم بتعبئة الجدول بالأسفل بالإضافة الى تقديم عرض تقني خاص بالشركة

<u>Traction machine</u>	Yes/No	Comment
Manufacture	Manufacture	
Model	Model	
Origin	Origin	
Delivery time	Delivery time	
3 years warranty and 10 years availability for spare parts	Yes No	
Applications: rehabilitation.	Yes/No	
For lumbar and cervical.	Yus/No	
User friendly	Yes/No	
User defined rest , treat time , hold .	YestNo	
Real time digital treatment tracking	Yus/No	
Therapy: continuous, intermittent, harmonized, progressive.	YesiNo	
Traction force: 0-90 kg.	Yes/No	
Force time: / rest time: up to 98s.	Yes/No	
Treat time: 0 to 90 min.	Yes/No	
Safety: audible signal (patient interrupt switch activated ,end of treatment)	Yes/No	
With all accessories: cervical traction attachment, pelvic attachment, knee bolsters, Flexion Stools, traction cord with metal fitting, patient interrupt switch, cervical harness separated bar, include in main price.	Ves/No	
Optional price :		
Manufacture of Examination bed with 3 years warranty	Manufacture	
Price for bed	Price	
model	Model .	
Origin	Origin	
	Y.es/No	
Examination bed for retraction		
Examination bed for retraction Price include user and service manual	Yes/No	

Infrared with ultraviolet

Application: rehabilitation.
Origin: USA, EURO.
Infrared with ultraviolet.
Adjustable arm and highest.

• يرجى التكرم بتعبئة الجدول بالأسفل بالإضافة الى تقديم عرض تقني خاص بالشركة

Infrared ultraviolet	Yes/No	Comment
Manufacture	Manufacture	
Model	Vlodel	
Origin	Origin	
Delivery time	Delivery time	
Two years warranty and 10 years availability for spare parts	Yos/No	
Applications : rehabilitation	Yes/No.	
Infrared with ultraviolet.	Yes/No	
Adjustable arm and highest .	Yes/No	
Price include user and service manual	Yes/No	
Price include installation and training in site .	Yes/No	



Laser Therapy

Origin: USA, EURO.

Wide and clear display.

Continuous and pulse operation.

Pre-programmed treatment protocols.

Display of energy in unit at joules or joules/cm².

Real time feedback of energy delivered.

Probe included.

Power:220v.

يرجى التكرم بتعبئة الجدول بالأسفل بالإضافة الى تقديم عرض تقني خاص بالشركة

Laser Therapy	Yes/No	Comment
Manufacture	Manufacture	
Model	Model	
Origin	Origin	
Delivery time	Delivery time	
3 years warranty and 10 years availability for spare parts	Yes/No	
Applications :physiotherapy.	Yes/No	
Wide and clear display.	Yes/No.	
Continuous and pulse operation.	Yes/No.	
Pre-programmed treatment protocols.	Yes/No	
Display of energy in unit at joules or joules/cm ² .	YestNo	
Real time feedback of energy delivered.	YestNo	
Probe included.	Yes/No	
Power:220v.	Ves/No	
Price include user and service manual	VesiNo	
Price include installation and training in site .	Yes/No	



Electrical stimulation Therapy

Application: : rehabilitation.

Origin: USA, EURO.

System include:

Electrical stimulation with all stander accessories to be ready to use

2 channel (independent)

Display: wide and good resolution screen.

User friendly and simple.

With standred accessories and electrodes.

• Kindly your offer will not be taken into our consideration if you don't fill below table.

Electrical stimulation Therapy	Yes/No	Comment
Manufacture	Manufacture	
Model	Model	
Origin	Origin	
Delivery time	Delivery time	
3years warranty and 10 years availability for spare parts	Yes/No	
The system include: rehabilitation.	Yes/No	
Electrical stimulation with all stander accessories to be ready to use	Yes/No	
2 channel (independent)	Vos/No	
Display: wide and good resolution screen.	Yes/No	
User friendly and simple.	Yes/No	
With standred accessories and electrodes.	Yes/No	
Price include user and service manual	Reice	
Price include installation and training in site .		
Any update of Software will be free of change if needed after of end the warranty.		



Ultrasound Therapy

Application: : rehabilitation.

Origin: USA, EURO.

System include:

Ultrasound with all stander accessories to be ready to use

Display: wide and good resolution screen.

Pulse and continues therapy .Variable output: from 16 to 100 Hz.

Head warmer for probe. User friendly and simple.

Different clinical protocol.

1 & 3 MHz ultrasound frequencies pulse and continuous mode. Sound head coupling features with both visual & audio feedback. Head of ultrasound probe $5cm^{-2}$.

Option price: Cart . Probe: 5cm².

* Kindly your offer will not be taken into our consideration if you don't fill below table.

Ultrasound	Yes/No	Comment
Manufacture	Manubieture	
Model	Model	
Origin	Origin	
Delivery time	Delivery time	
3 years warranty and 10 years availability for spare parts	Yes/No	
Applications :	Yes/No	
The system include: rehabilitation.	Mes/No	
Ultrasound with all stander accessories to be ready to use	YescNo	
Display: wide and good resolution screen.	YeseNo	
Pulse and continues therapy .	Yes/No	
Variable output : from 16 to 100 Hz.	Yes/No	
Different clinical protocol.	Yes/No	
1 & 3 MHz ultrasound frequencies pulse and continuous mode.	Yes/No	
Sound head coupling features with both visual & audio feedback.	Yes/No	
Head of ultrasound probe 5cm ² . Included in u.p	Yesthan	
Option price:		
Cart .	Price	
Probe: 5cm ² .	Price	
Price include user and service manual	Price	
Price include installation and training in site .		
Any update of Software will be free of change if needed after of end the warranty.		

